

-the Better Choice

Great White® Platinum Stirrups

*Gynecologic, laparoscopic
& urologic procedures
have never been easier*

The quick-connect, blade-mounting style with extended lateral boot system can be locked into any position within the entire lithotomy and abduction/adduction range

Features

- *Raise, lower, abduct, adduct quickly and easily with the Trigger-Grip Handle*
- *Extended lateral boot design helps reduce the risk of peroneal nerve injury*
- *600 lb. surgical patient weight capacity*
- *Deluxe boot liners wrap around the foot, ankle and calf providing safe, comfortable and secure patient positioning*
- *Lithotomy Range -35° low to 90° high; Adduction -9° and abduction +25°*
- *Length and lithotomy indicators help to position and reposition the legs*
- *Rotates on three axes*

#800-0007 Great White® Platinum Stirrups

Required Accessories

#800-0127 Spring Loc Clamp (2 required)

Optional Accessories

508-0113 Great White® Platinum Boot Pads

#800-0160 Clear View Leg Drapes

#800-0075 Dolly, Great White® Stirrup

#800-0284 Stirrup Wall Rack



*Increase your
surgeons' O.R.
performance*



Great White  Platinum Stirrups

- 600 lb (272 kg) surgical patient weight capacity
- Boot free-floats throughout each range reducing the risk of superficial nerve injury
- Raise, lower, abduct, adduct quickly and easily with the Trigger Grip handle
- Lithotomy and length indicators help to position and reposition the legs
- Deluxe boot stirrup liners wrap around the foot, ankle and calf providing safe, comfortable and secure patient positioning

Increase your surgeons' operating room performance with the Great White Platinum Stirrups which permits simple adjustment of lithotomy and abduction intraoperatively. A simple squeeze of the Trigger Grip trigger provides complete control of lithotomy and abduction/adduction without compromising the sterile field. Simply release the trigger on the handle to securely lock the stirrup in all directions. Quick-connect, blade-mounting style with extended lateral boot system, the Great White Platinum Boot Stirrups' range of motion is -35° low lithotomy to 84° high lithotomy. Adduction -9° and abduction 25°, with the control to lock the stirrup anywhere in between. The Great White ships fully assembled. Boot liners are included. Two required SpringLoc Clamps sold separately.

#800-0007 Great White Platinum Stirrups

Required Accessory: Purchase Two

#800-0127 SpringLoc Clamp

Optional Accessories


#800-0075 Dolly, Great White Stirrup

#800-0284 Stirrup Wall Rack



Organize your accessories with the Stirrup Dolly (left) and Stirrup Wall Rack (right).

Optional Accessories



Attaches to any general surgical table to contain fluids for hysteroscopy or cystoscopy procedures.

#800-0096 Uro-Trapper System



Allows visibility of patient's leg during intra-operatively adjustments. 10/cs.

#800-0160 Disposable Clear View Leg Drapes, Sterile



Prevent patients from sliding during steep reverse Trendelenburg.

#800-0055 Nissen Straps