



## Great White **Maxima** Stirrups

- 800 lb (363 kg) surgical patient weight capacity
- Boot free-floats throughout each range reducing the risk of superficial nerve injury
- Raise, lower, abduct, adduct quickly and easily with the Trigger Grip Handle
- Lithotomy and length indicators help to position and reposition the legs
- Extended lateral boot design helps reduce the risk of peroneal nerve injury

Increase your surgeons' operating room performance with the Great White Maxima Stirrups which permit simple adjustment of lithotomy and abduction intraoperatively. A simple squeeze of the Trigger Grip Handle provides complete control of lithomy and abduction/adduction without compromising the sterile field. Quick connect blade-mounting style with extended lateral boot system. Deluxe boot stirrup liners wrap around the foot, ankle and calf provid-

ing safe, comfortable and secure patient positioning. Great White Maxima Stirrups range of motion is -35° low lithotomy to 84° high lithotomy. Adduction -9° and abduction 25°, with the control to lock the stirrup anywhere in between. The Great White ships fully assembled. Boot liners are included. Two required SpringLoc Clamps sold separately.

**#800-0245** Great White Maxima Stirrups

*Required Accessory: Purchase Two*  
**#800-0127** SpringLoc Clamp

*Optional Accessory*

**#800-0075** Dolly, Great White Stirrup