

Schure Loc™ XPS

Easy as 1, 2, 3

- *Squeeze and hold lever to release all device joints*
- *Position limb in exact position*
- *Release lever to securely lock into position*

The Schure Loc™ XPS combines strength and flexibility to provide optimal intraoperative positioning for all open and arthroscopic limb procedures.

State-of-the-art patient positioning device improves OR staff performance

Schure Loc™ XPS

- *Securely and safely holds patient arm in any desired position during arthroscopic surgery*
- *Compact and lightweight*
- *No foot pedals, hoses, hydraulics or power cords*
- *Ships assembled and ready for immediate use*
- *Single lever releases all joints for positioning and relocks when released*
- *Decreases surgeon fatigue and improves intraoperative control of surgical limb*
- *Schure Loc™ XPS reduces staffing needs providing a more cost-effective procedure*
- *Schure Loc™ XPS Shoulder Kit conveniently includes a sterile drape and foam hand/forearm holder, and self-adherent wrap*

Schure Loc™ XPS currently outperforms similar equipment for shoulder surgery in both the beach chair and lateral position. The Schure Loc™ XPS is the perfect companion piece for any beach chair positioner—easily attaching to any OR table rail while freeing your surgical assistant to do much more than steady the patient's limb. Get the most from your team by adding this high-quality, easy-to-use device to your shoulder surgery equipment.

#800-0280 Schure Loc™ XPS

Optional Accessory

#800-0281 Disposable Schure Loc™ XPS Sterile Kit, 10/cs

Includes Sterile Forearm Wrap, Sterile Self-Adherent Coban™ and Sterile Drape.



Provides excellent positioning and surgical access for beach chair, lateral shoulder and supine elbow procedures.